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
#### FLEXIBLE THERMOMETER

For adults and grown up children:

- for rectal use (in the anus),
- for oral use (in the mouth),
- for axillary use (under the arm).


- Dimensions: 133mm (L) x 21mm (l) x 12mm (H)
- Weight: Approx. 9.5g (battery included)
- Battery : Buttoncell 1.5V (LR41)

#### SOOTHER THERMOMETER

Especially for babies:  
 Oral use (in the mouth) / Temperature in 2-4 min.  
 Wait 6-8 min between each measurement.

- Dimensions: 40mm (L) x 43mm (l) x 53mm (H)
- Weight: Approx. 14g (battery included)
- Teat size : Medium
- Battery: DC 1.55/1.5v (SR41/LR41)
- Power consumption: 0.2mW

#### BATTERY REPLACEMENT

- Remove the battery cover by lifting it up with a pointed non-conductive object.
- Gently pull out the circuit board with the battery compartment about 1cm.
- Use a pointed non-conductive object again to dislodge the battery by pushing it out of its location. Dispose of the used battery in accordance with the standards and the law in force.
- Insert a new buttoncell battery (LR41) with the positive (+) side facing up.
- Reinsert the battery compartment into position and replace the battery cover.

#### BATTERY REPLACEMENT

The battery cannot be extended for the soother. Dispose of the thermometer responsibly. Do not dispose it with household waste, bring it back to a collection point for e-waste.

#### SPECIFICATIONS

- Display unit: 0.1°C
- Temperatures are measured in degrees Celsius (°C). The symbol is displayed on the upper right corner of the screen.
- Measuring range: 32.0°C-42.9°C
- Accuracy: ±0.1°C between 35.5°C-42°C and 18°C-28°C for ambient temperature range. ±0.2°C for other measurements and operating temperature range.
- Operational mode: Direct mode
- Display: 3.5 digit LCD display
- Memory: Save the last measured value in memory. CAUTION: It is displayed momentarily after switching on the thermometer and is accompanied by the symbol **M**.
- Expected lifetime: More than 200 hours of operation. (Approximately 2 years if used 10 minutes per day). If you don't manually turn off the thermometer, it will turn off automatically after a few minutes.
- Operating ambient temperature range: Temperature: 5°C - 40°C / Relative humidity: 15% - 95%RH / Atmospheric pressure: 700hPa - 1060hPa
- Classification : Type BF
- Ingress protection rating: IP22

#### TROUBLESHOOTING

The displayed temperature is below 32.0°C. Turn off the thermometer. Wait a minute and switch on again to take a new measurement.

The displayed temperature is higher than 42.9°C. Turn off the thermometer. Wait a minute and switch on again to take a new measurement.

The system does not work properly. Remove the battery. Wait a minute, replace the battery and then turn the thermometer back on.

When the battery is too low, the symbols display flashing. Replace the battery.

### IMPORTANT! READ CAREFULLY AND KEEP FOR FUTURE REFERENCE.

These digital thermometers provide a quick and extremely accurate reading of body temperature. To better understand their operation and to enjoy many years of reliable use, please read all instructions before using it.

#### WARNINGS!

- Choking Hazard: The flexible thermometer's cap and battery may be fatal if swallowed.
- Do not allow children to use these devices without adult supervision. Store them out of the reach of children.
- Do not use the flexible thermometer in the ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
- Keep away from electrostatic or electromagnetic waves (ex: TV, microwave or telephone).
- Do not expose the thermometers to direct sunlight, high temperature, humidity or dust.
- Do not expose the battery to a high heat source as it may explode.
- Note: Use of the cover when the thermometer is wet. This may result in inaccurate readings.
- It is recommended to check the performances every two years. Stop using if broken.
- Remove battery from the flexible thermometer when not in operation for a long time.
- The use of temperature readings for self-diagnosis is dangerous. Consult your doctor for the interpretation of results. Self-diagnosis may lead to the worsening of existing disease symptoms.
- Do not attempt to take the temperature when the thermometer is wet. This may result in inaccurate readings.
- Do not bite, bend or drop. This could break the thermometers and/or cause injury.
- Do not attempt to disassemble or repair the thermometers. This may result in inaccurate readings.
- After each use, disinfect the thermometers especially in case of the device is used by more than one person.
- Do not force the flexible thermometer into the rectum. In case of pain, remove the thermometer and stop taking the temperature. Follow this instruction to avoid risk of injury.
- Do not use the flexible thermometer orally (for baby use) rectally.
- For children who are two years old or younger, please do not use the devices orally.
- When the thermometer has been exposed/stored at a temperature below 5°C (41°F) and above 40°C (104°F), leave it at room temperature for about 15 minutes before using it.
- Inspect carefully before each use. Pull the soother in all directions. Throw away at the first signs of damage or weakness.
- Never attach other ribbons or cords to a soother, your child may be strangled by them. Stopping it at the slightest sign of deterioration or fragility.
- Do not use the soother thermometer as a regular soother.

#### CLEANING AND DISINFECTION

Follow cleaning instructions. Do not wash in any way other than those indicated. Always keep the probe and teat clean. Use a damp cloth with warm water and a little mild soap. Rinse carefully without immersing them (thermometers are not waterproof), then wipe with a soft, dry cloth.

- To disinfect the probe, use an agent based on isopropyl alcohol or 75% ethanol.
- Never steam sterilize or boil the thermometer, do not wash in the dishwasher, do not put in a microwave sterilizer or in the microwave.
- Do not clean the thermometers with benzene-based solvent, thinner, gasoline, or other aggressive agent.
- Do not immerse the thermometer probe in alcohol or hot water (> 50°C)
- Do not use an ultrasonic system to clean the thermometers.

#### WARRANTY

2 years warranty on presentation of the purchase invoice.  
 Does not cover battery and thermometer damage caused by improper use.  
 Contact after-sales service on 01 49 15 96 40.

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#### TEMPERATURE

The human body temperature is normally between 36.5 and 37.5°C. This temperature can vary depending on age, activity and time of day, etc. Also, temperature varies depending on how it is measured\*. When you report your temperature to a professional, it is advisable to indicate the method of measurement.


\* A rectal measurement gives the most accurate temperature. Oral measurements are approx. 0.5°C lower than the actual body temperature and the temperature under the arm is 0.5°C lower than in the mouth. In the morning, your temperature is lower than in the afternoon.

#### TEMPERATURE MEASUREMENT :

- Press the On/Off button located next to the LCD display. A tone will sound and the screen shows 88.8°C followed by the last recorded temperature. After showing the self-test temperature, the thermometer is ready to take the temperature.
- Depending on the thermometer used, position it according to the place chosen:
  - mouth, rectum or armpit for the flexible thermometer.
  - mouth for the soother thermometer.

#### A. Oral temperature:

Flexible thermometer:  
 Place the thermometer under the tongue in the position shown (Fig 1). Close your mouth and breathe normally and regularly through your nose to prevent the measurement from being influenced by the passage of exhaled/inhaled air.

#### Soother thermometer:

Place the soother thermometer and its silicone teat in the mouth (Fig 2). Close your mouth and breathe normally and regularly through your nose to prevent the measurement from being influenced by the passage of exhaled/inhaled air.  
 Note: Drinking or eating before can distort the reading, wait 30 minutes before taking the temperature.

#### B. Rectal temperature measurement:

Lubricate the silver end of the probe with petroleum jelly to facilitate insertion. Gently insert the thermometer probe into the rectum for a length of 1cm.

#### C. Axillary temperature measurement:

Wipe the armpit dry. Place the probe under armpit and keep your arm pressed firmly at the side of the body. From a medical point of view, this method always gives inaccurate readings. It should therefore not be used if precise temperatures are required.

- The degree °C symbol flashes throughout the temperature measurement. When it stops flashing, an alarm sounds for approximately 10 seconds. At the same time, the measured temperature is displayed on the LCD display. The minimum duration from the temperature measurement until the alarm (beep) must be respected without exception. The measurement continues even after the alarm signaling the end of the measurement. To obtain a better measurement of body temperature, we recommend leaving the probe in the mouth or rectum for a minimum of 2 minutes and 5 minutes under the armpit, ignoring the alarm. A minimum measurement interval of 30 seconds must be observed.
- To extend battery life, press the On/Off button to turn off the thermometer when finished taking the temperature. If you do not press this button the thermometer will automatically turn off after 8 minutes.

\* Note: Normally the alarm sounds "Bi---Bi---Bi" and accelerates when the temperature reaches 37,8°C or more. The alarm then sounds "Bi-Bi-Bi---Bi-Bi-Bi".

Acute Idecs Co. Ltd.  
 3F No. 11 Lane 35, Jihu Road, Neihu District Taipei 114, Taiwan  
 Ma de in China

WellKang Ltd., The Enterprise Hub, NW Business Complex, 1 Beraghmore Rd. Derry, BT48 8SE, Northern Ireland. Tels: +353 (0)144 33 560

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